

LOW-VACUUM WOUND DRAINAGE

Closed drainage system with manual low-vacuum regulation for a gentle drainage of blood and secretion following e.g.

- dissection of axillary lymph nodes
- mammoplasty

FIELDS OF APPLICATION:

- Abdominal Surgery
- Mamma Surgery
- Plastic Surgery
- Neurosurgery



SECURO DRAIN® silicone reservoir for manual vacuum regulation

FEATURES:

- the non-adhesive silicone material allows an easy, pain-avoiding drainage removal
- the ribbed inner wall design guarantees a problem-free drainage, even when the drainage tube is collapsed
- internal anti-reflux valve prevents backflow of fluid to the patient
- transparent reservoir and graduations make it easy to identify and measure exudate

Description	Ref.No.	Length of Perforation	Size	Packaging Unit	Drainage-Set with 150cc-Reservoir	Packing Unit
Type: flat, 3/4 perforation	01.3307	15 cm (3/4)	7 mm	20 EA/Box	01.3357	5 EA/Box
	01.3310	15 cm (3/4)	10 mm	20 EA/Box	01.3360	5 EA/Box
Type: flat, complete perforation	01.3404	15 cm (4/4)	4 mm	20 EA/Box	01.3454	5 EA/Box
	01.3407	20 cm (4/4)	7 mm	20 EA/Box	01.3457	5 EA/Box
	01.3410	20 cm (4/4)	10 mm	20 EA/Box	01.3460	5 EA/Box

DRAIN , round

Description	Ref.No.	Length of Perforation	Total Length	Packaging Unit
Type: round	<i>10FR</i>	23 cm	125 cm	20 EA/Box
	<i>15FR</i>	23 cm	125 cm	20 EA/Box
	<i>19FR</i>	23 cm	125 cm	20 EA/Box

Reservoir

Description	Ref.No.	Contents	Packaging Unit
Reservoir	01.3300	150 cc	5 EA/Box
	01.3303	400 cc	2 EA/Box

STERILE EO LATEX-FREE 0482

SPECIFICATIONS:

- Material:
Drain:..... Silicone, radiopaque
- Dimensions Flat Drain:..... 80cm (total)
Flat part:..... 20cm
Extension line:..... 60cm