



SUCTION IRRIGATION SYSTEM

SUCTION IRRIGATION SYSTEM



SUCTION IRRIGATION SYSTEM

- Color-coded buttons and tubings for easy identification
- Various probes meet all kinds of demands
(available in 28cm, 33cm, 45cm working length and 5mm, 10mm diameter)

- One way valve prevents the saline from refluxing
- Ergonomically designed handpiece can be operated by either hand
- Extra smooth movement of the handpiece prevents fatigue
- Clear handpiece for unobstructed visualization
- Quick disconnection allows easy exchange of instrumentation
- Chamfered probe design reduces damage

ELECTRODES WITH SUCTION CHANNEL



SUCTION IRRIGATION SYSTEM

DESCRIPTION	PACKAGING
Suction irrigation set 5mm/33cm	10 sets/box
Suction irrigation set 5mm/28cm	10 sets/box
Suction irrigation set 5mm/45cm	10 sets/box
Suction irrigation set 5mm/33cm with L hook	10 sets/box
Suction irrigation set 5mm/33cm with spatula	10 sets/box
Suction irrigation set 5mm/33cm with button electrode	10 sets/box
Disposable probe for S/I set 5mm/28cm	20 pcs/box
Disposable probe for S/I set 5mm/33cm	
Disposable probe for S/I set 5mm/45cm	
Disposable probe for S/I set 10mm/33cm	
Reusable probe for S/I set 5mm/33cm	20 pcs/box
Reusable probe for S/I set 5mm/45cm	
Reusable probe for S/I set 10mm/33cm	
Monopolar L hook for S/I set 5mm/33cm	
Monopolar spatula for S/I set 5mm/33cm	
Monopolar button electrode for S/I set 5mm/33cm	



**ENDO
POUCH SERIES**

ENDO POUCH SERIES



ENDO POUCH SERIES



Endo pouch with introducer



Endo pouch with suture tie



Endo pouch with memory wire

BENEFITS OF DETACHABLE ENDO POUCH

Unimax detachable endo pouch is a patented device which provides a quicker and easier application of endo pouch. This unique design allows the pouch to be detached from the introducer and left within Patient body during procedure. Consequently, it increases the convenience of procedure and shortens the procedure time.



- Unique patented design ensures ease of use
- 10mm introducer provides ease of insertion
- Self-opening pouch design provides easy and safe specimen capture
- High strength PU pouch material assures no leakage
- Enhanced metal ring ensures large pouch opening and keeps the pouch in place when loading specimen

ENDO POUCH SERIES

DESCRIPTION	PACKAGING
Endo pouch with introducer 3"*6"	20 pcs/box
Endo pouch with introducer 5"*7"	20 pcs/box
Endo pouch with suture tie 3"*6"	20 pcs/box
Endo pouch with suture tie 5"*7"	20 pcs/box
Endo pouch with memory wire 3"*6"	10 pcs/box
Endo pouch with memory wire 5"*7"	10 pcs/box
Detachable endo pouch 3"*6"	5 pcs/box
Detachable endo pouch 5"*7"	5 pcs/box



MONOPOLAR SCISSORS


MONOPOLAR SCISSORS



MONOPOLAR SCISSORS

- Stainless steel blades sustain sharpness throughout the procedure
- Precise & sharp blades for smooth cuttings

360° rotation wheel allows the shaft to rotate for 360 degrees



Extended PTFE insulation reduces the chance of thermal damage

- Ergonomic handle increases comfort during operation
- Extra smooth manipulation

MONOPOLAR SCISSORS

DESCRIPTION

Monopolar scissors 5mm/33cm


PACKAGING

5 pcs/box



LAPAROSCOPIC
GRASPING
INSTRUMENTS

LAPAROSCOPIC GRASPING INSTRUMENTS



Duckbill grasper



Grasper



Clinch



Rat tooth



Dissector



Babcock

GRASPING FORCEPS

PRODUCT DESCRIPTION

PACKAGING

Straight grasping forceps 5mm/33cm

5 pcs/box

Duckbill grasping forceps 5mm/33cm

5 pcs/ box

Maryland dissecting forceps 5mm/33cm

5 pcs/box

Clinch grasping forceps 5mm/33cm

5 pcs/box

Babcock grasping forceps 5mm/33cm

5 pcs/box

Rat tooth grasping forceps 5mm/33cm

5 pcs/box