

MEDI-I-TAPE®

INSTRUMENT IDENTIFICATION TAPE

MEDI-I-TAPE® color coding is a fast, easy and safe way to identify surgical instruments

MEDI-I-TAPE®

- allows quick identification and sorting of instruments
- marks instruments used only in special surgical applications or surgical units
- helps to assign instruments to the different surgical kits



9 colors allow numberless color coding combinations

FEATURES:

- can be washed and disinfected
- can be repeatedly autoclaved or sterilized by ETO
- extra long (7.6m / roll)
- supplied in a reclosable plastic box

Description	Ref.No.	Length	Width	Color
MEDI-I-TAPE®	01.4001	760 cm	6.5 mm	White
	01.4002	760 cm	6.5 mm	Blue
	01.4003	760 cm	6.5 mm	Green
	01.4004	760 cm	6.5 mm	Red
	01.4005	760 cm	6.5 mm	Yellow
	01.4006	760 cm	6.5 mm	Purple
	01.4007	760 cm	6.5 mm	Orange
	01.4008	760 cm	6.5 mm	Brown
	01.4009	760 cm	6.5 mm	Black

LATEX-FREE A CE

SPECIFICATIONS:

- Material:
Tape:PVC / white rubber resin
- Packaging:.....single packed, non sterile
Unit:.....1 roll / box