

Special Oil Spray

FOR THE MAINTENANCE OF
SURGICAL INSTRUMENTS

Special Oil Spray is used for the maintenance of surgical instruments with hinges, like scissors and forceps. Due to the employment of paraffinum liquidum the product is toxicologically harmless and offers the highest degree of purity. Special Oil Spray preserves the value of instruments by providing an optimal lubrication and an effective corrosion inhibition.



FEATURES:

- prolongs the shelf life of surgical instruments
- avoids incrustation on instrument surfaces
- reduces protein- and mineral deposits on the instrument surface
- minimizes stains and spots

Description	Ref.No.	Contents
Special Oil Spray	01.6055	500 cc

LATEX-FREE 

SPECIFICATIONS:

- Material:
Ingredients: Paraffinum liquidum
Propellant: Butane / Propane (CFC-free)
- Packaging:
Unit: 12 cans / box

Silicone Spray

FOR THE MAINTENANCE OF
RUBBER, LATEX AND PLASTIC
PRODUCTS

For the maintenance of medical instruments and all kinds of steel goods including rigid endoscopes, turbines and hinges. Silicone Spray is toxicologically harmless, thus being an excellent lubricant for flexible endoscopes, catheters and drainages.



FEATURES:

- guarantees to an optimal lubrication of mobile instrument parts
- helps to protect areas of mucous membrane during insertion of drainages or indwelling catheters
- offers an effective protection against corrosion
- prevents thermolabile materials like latex or rubber from swelling
- suitable for the protection of endoscope valves
- **WARNING:** Silicone Spray should not be used for silicone products (suitable: SPECIAL OIL SPRAY)

Description	Ref.No.	Contents
Silicone Spray	01.3055	500 cc

LATEX-FREE 

SPECIFICATIONS:

- **Material:**
Ingredients: Dimethylpolysiloxan
Propellant: Butane/Propane (CFC-free)
- **Packaging:**
Unit: 12 cans / box