




Introduction to Inara





Advantages of Inara = II

- Inara comes ready for insertion in insertion tube, therefore less number of steps while insertion than other IUDs
- Nylon Suture is softer and hence is more comfortable for the couple
- Superior Film Film pouch keeps the copper parts bright
- Geometric design of Inara's frame ensures lower expulsion rate



Advantages of Inara = III

- Geometrical structure and insertion method lowers the probability of accidental perforation of fundus
- Printed insertion tube adds convenience and accuracy for providers
- Inara range of IUDs offer long term protection (5 years for Cu 375 & Cu 375 Sleek and 3 years for Cu 250)




Additional Advantages of Inara Cu 375 Sleek

- Shorter length of Inara Cu 375 Sleek allows usage by women with smaller uterus (particularly in Asian region). Small length also provides more comfort for user. (Smaller by 6.2mm than Cu 375)




Additional Advantages of Inara Cu 250

- Lower copper content in Inara Cu 250 allows use by women with anemic or excess bleeding condition.
- Inara Cu 250 allows spacing of up to 3 years




Introduction to Pregna TCU 380A




Advantages of Pregna TCu 380A Part - I

- Special modification in T-frame to increase its strength.
- It is packed in Film-Film pouch. Film-Film pouch is better in protecting Copper parts of TCu 380A and is far superior in moisture resistance than Tyvek Film Pouch. It also makes the whole packaging look neater and more professional.
- Solid rod of Pregna's TCu 380A has a special fin design, which holds the solid rod in place even while the service provider is holding the loaded insertion tube vertically.



Advantages of Pregna TCU 380A Part - II

- Nylon suture is softer and hence is more comfortable for the couple.
- Printed insertion tube adds convenience and accuracy for the provider
- Pregna TCU 380A offers long term protection of up to 10 years
- Being a non-barrier method, offers complete intimacy for both the partners
- Fertility is restored on removal and absolutely safe for lactating mothers



Advantages of Pregna TCU 380A Part - III

- Provides complete freedom from consistent follow up & motivation
- High efficacy of more than 99%
- Insertion & removal of Pregna TCU 380A is simple and safe

Introduction to Safeload

Unique Safeload by Pregna




Advantages of Safeload - II

- The loading of 'T' frame is accurate and requires minimum effort
- Both the arms of T-frame are evenly loaded, allowing correct insertion and eliminates chances of rupturing the pouch.
- Offers superior no touch technique
- Ensures that the service provider's fingers remain agile and fatigue free
- Printed insertion tube adds convenience and accuracy for the provider

Advantages of Safeload - II


- Nylon suture is softer and hence is more comfortable for the couple
- Modified insertion tube ensures fundal placement
- Superior Film Film pouch ensures bright copper parts
- Pre-loaded solid rod reduces one step in the insertion method, hence more convenience for provider
- Safeload also has all the advantages of Pregna TCu 380A



Introduction to Silverline

Advantages of Silverline - I


- Silverline's flexible arms offer a perfect fundal placement
- Silverline is gentle during insertion and removal




Advantages of Silverline - III

- In Silverline the copper wire has Silver Core, which prevents copper wire from fragmenting. Dissolution of wire stops once it reaches the Silver Core, as silver does not dissolve under contact with uterine fluids

Silverline 380Ag



Silverline 200Ag



Advantages of Silverline - III

- Silverline offers a long term protection (5 years for 380Ag and 3 years for 200Ag)
- Silverline requires 1/3rd dilation of cervix during insertion, in comparison to Copper T380A
- Nylon Suture is softer and hence is more comfortable for the couple
- Superior Film Film pouch keeps the copper parts bright and shiny
- Printed insertion tube adds accuracy and convenience for providers

Additional Advantages of Silverline

200Ag

- Lower copper content in Silverline 200Ag allows use by women with anemic or excess bleeding condition
- Silverline 200Ag allows spacing of up to 3 years

Silverline
Intrauterine Contraceptive Device

Thank You

