

Silicon Felt

SCIBA

They are designed for use in several and special-invasive surgery cases i.e. tympanoplasty, septoplasty, cardiosurgery e.t.c. and are available at several sizes, thicknesses and hardness levels.

They are available in sterile or non sterile individual packaging. They are silicon of high purity and quality level. They can stay in the body for up to 30 days.



TYPE	HARDNESS	THICKNESS	DIMENSION
SILICON FELT S4	37 shore	0.12mm-0.25mm	50 X 60mm
		0.12mm-0.25mm	70X90mm
		1.00mm-1.50mm	70X180mm
			90X150mm
			150X160mm
			180X180mm
SILICON FELT S6	52 shore	0.12mm-0.25mm	50 X 60mm
		0.12mm-0.25mm	70X90mm
		1.00mm-1.50mm	70X180mm
			90X150mm
			150X160mm
			180X180mm
SILICON FELT S7	66 shore	0.12mm-0.25mm	50 X 60mm
		0.12mm-0.25mm	70X90mm
		1.00mm-1.50mm	70X180mm
			90X150mm
			150X160mm
			180X180mm
SILICON FELT S9 (reinforced)	90shore	0.12mm-0.25mm	50 X 60mm
		0.12mm-0.25mm	70X90mm
		1.00mm-1.50mm	70X180mm
			90X150mm
			150X160mm
			180X180mm

Ordering Information :

Felt dimension upon request, please contact us for complete list of references

Our product is sterile, single use packaged, with no requirements for special storing conditions

180X180mm

CE