

# Aortic Punch

**SCIBA SA**

Aortic Punch is a sterile, single-patient use disposable aortic punch designed to produce precise dimensional circular or oval openings in cardiac surgery during coronary artery bypass crafting procedures. The product is available in 16 different sizes to meet any surgeon's needs



Code #	Punch Size	Assembly
APL28	2,8 mm	Long
APM28	2,8 mm	Medium
APL36	3,6 mm	Long
APM36	3,6 mm	Medium
APL40	4,0 mm	Long
APM40	4,0 mm	Medium
APL44	4,4 mm	Long
APM44	4,4 mm	Medium
APL48	4,8 mm	Long
APM48	4,8 mm	Medium
APL52	5,2 mm	Long
APM52	5,2 mm	Medium
APL56	5,6 mm	Long
APM56	5,6 mm	Medium

## Ordering Information :

Felt dimension upon request, please contact us for complete list of references

Our product is sterile, single use packaged, with no requirements for special storing conditions

Special packaging /volume upon request

**CE**