



Features

- Stainless steel semi sharp cutting edge for soft or weak tissue
- Size embossed on each curette for easy identification
- Wide range of sizes
- Sterile and ready to use
- Disposable

Used in the following minor skin surgical procedures






- Curettage of a basal cell carcinoma
- Curettage of a verruca
- Curettage of other benign neoplasms of the skin

Clinical areas

- Dermatology
- Surgery
- Podiatry
- Veterinary

Package

- 20 curettes in a dispenser box
- 20 boxes per shipping carton

REF	SIZE	
MK402	2 mm	
MK403	3 mm	
MK404	4 mm	
MK405	5 mm	
MK407	7 mm	

Actual size



A dermal curette biopsy technique for various dermatological diseases



A punch (Trepan) biopsy technique for various dermatological diseases

Product specifications may be changed without notice

