

# ELECTROSURGERY

## ELECTROCAUTERIES

FIAB bipolar sterile cautery pens feature perfect haemostasis with fast and effective precision cauterization of small blood vessels. Being powered by an internal long life battery, they are extremely easy to handle; several models are available, different from each other for working temperatures (1200°C, 800°C and 600°C), shape, thickness and length of the cautery tip. The available range allows to choose a cautery pen that perfectly matches the clinical demands of the procedure.



### Disposable electrocauterics

#### High temperature


F7244	28 mm fine tip
F7266	28 mm large tip
F7234	75 mm fine tip
F7277	125 mm fine tip

#### Medium temperature

F7288	28 mm thick tip
-------	-----------------

#### Low temperature

F7255	28 mm fine tip
F7255UF	28 mm ultra fine tip



### DISPOSABLE ELECTROCAUTERIES - QUICK REFERENCE GUIDE

MODEL	F7234	F7244	F7255	F7255UF	F7266	F7277	F7288
TIP SHAPE	FINE	FINE	FINE	ULTRA FINE	LARGE	FINE	THICK
TIP LENGTH (mm / inch)	75 / 3,0	28 / 1,1	28 / 1,1	28 / 1,1	28 / 1,1	125 / 5,0	28 / 1,1
TEMPERATURE AT THE ACTIVATION (°C / °F)	1200 / 2200	1200 / 2200	600 / 1100	600 / 1100	1200 / 2200	1200 / 2200	800 / 1400
ACTIVATION PRESSURE (gr)	200	200	200	200	200	200	200
POWER SUPPLY (V)	3,0	3,0	1,5	1,5	3,0	3,0	3,0
DIAMETER (mm)	16	16	16	16	16	16	16
WEIGHT (gr)	78	65	45	45	65	78	65


## *Reusable electrocauteries*

### Reusable electrocauteries

- F7297 low temperature (600°C/1100°F)
- F7299 high temperature (1200°C/2100°F)

### Disposable sterile tips

- F7298F fine tip, 30 mm, suitable for F7297 and F7299
- F7298L large tip, 30 mm, suitable for F7299
- F7298FG long fine tip, 120 mm, suitable for F7299



## NERVE STIMULATOR

### Neuropacer

- F1745 Disposable, sterile, battery operating, complete with active electrode and grounding needle

