

Sternum Closure System

for sternum closure after
intrathoracic procedures



packaging for single pack and multi-pack



Sternum wire, round packed

FEATURES:

- optimal geometric needle cut facilitates an easy penetration of the sternum
- optimized alloy for highest tensile loading
- implant quality

Description	Ref.No.	Length	Wire Ø	metric / USP	Wires / Foil	Packaging Unit
Sternum Closure System with Lancet-Trocar-Needle HSLX 50, round packed	49.F71R	75 cm	0.78 mm	7 / 5	1 EA	12 Foils / Box
Sternum Closure System with Lancet-Trocar-Needle HSLX 50, long packed	49.F71	75 cm	0.78 mm	7 / 5	1 EA	24 Foils / Box
	249.F71	75 cm	0.78 mm	7 / 5	2 EA	24 Foils / Box
	449.F7	50 cm	0.78 mm	7 / 5	4 EA	24 Foils / Box
	249.F81	75 cm	0.85 mm	8 / 6	2 EA	24 Foils / Box
	449.F8	50 cm	0.85 mm	8 / 6	4 EA	24 Foils / Box
	449.F9	50 cm	0.95 mm	9 / 7	4 EA	24 Foils / Box

STERILE | EO | LATEX-FREE | 0482

SPECIFICATIONS:

- Material:
 - Wire:..... Stainless Steel, monofilament, implant quality
 - Needle:..... Stainless Steel
- Dimensions:
 - Needle:..... 50 mm curve length
 - Type:..... 1/2 circle, cutting, lancet-trocar-point cut
- Packaging:..... double packed, sterile